

# EOSedge™, A unique X-Ray imaging modality system



Ultra-low dose  
radiation



Versatility across  
imaging suite



High resolution  
images



High patient  
throughput




# EOSedge Referrers Guide



EOSedge™ System provides low dose, full body, stereo radiographic images of the patient in a functional position

**Exam Time:**  
4 Min



Stereos Workstation enables the conception of patient-specific models and reports based on EOSedge images

**Processing Time:**  
Basic reports 5-10 Min  
Advanced Reports 10-30 Min



## Clinical Applications

### Spine

- **Paediatric:** Scoliosis, Postural disorders
- **Adult:** Postural disorders, deformative and degenerative spine pathologies

### Pelvis and lower limbs

- Pelvis deformities and postural disorders
- **Lower limb deformities:** Leg length discrepancies and alignment disorder
- Management of hip and knee arthroplasties

## Typical Orders

Full Body

- Postural assessment
- Lower limbs alignment

Lower Limbs

- Lower limb alignment
- Post-op THA

Pelvis + Lumbar spine sitting

- Pelvic assessment

Spine  Brace  Brace off

- Postural Assessment

Spine Sitting

Spine Bending

Spine Flexion/extension